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drive/folders/1q-86qGK3xqhEvCPqCVJcvfRvQKc3t1-a

Pneumask PPE
Instructions for Use
Intended Use

This device is intended to be worn by health care professionals to protect the user from transfer of microorganisms, body fluids, and particulate material. It is intended to be a reusable substitute to an N95 respirator and face shield/eye goggle combination, or a PAPR helmet when no other appropriate means of protection exists due to resource limitations from pandemic conditions in the region of use.

Warnings

- ALWAYS remove masks very gently.
- NEVER pull directly away from your face, it can spring and damage your eyes.
- ONLY choose the filters according to the instructions in the “[Appropriate Filter Choice](#)” session.
- ONLY wear Pneumask PPE for protecting healthcare workers in a standard hospital setting. It is not appropriate for situations when protections against hazards other than bacteria and viruses are required.
- NEVER seal or block the chin valve.
- NEVER wear Pneumask PPE ALONE when a sterile field is required, unless the chin valve is sealed.
- If you are diagnosed to be claustrophobic, or have a history of panic disorder, using a full-face mask might trigger a panic attack.
- Individuals with a compromised respiratory system, such as asthma or emphysema, should consult a physician and must complete a medical evaluation prior to use.

Cautions

- Failure to follow all instructions and limitations on the use of Pneumask PPE and/or failure to wear Pneumask PPE during all times of exposure can reduce filtering effectiveness and may result in sickness or death.
- Leave the contaminated area immediately and contact your supervisor if dizziness, irritation, or other distress occurs.
- Leave the contaminated area immediately if the eye chamber becomes fogged or covered in condensation.
- Conduct a fit test according to the instructions in the “[Fit Test Instructions](#)” section. If you cannot pass the fit test, do not use the mask.
- Conduct a user seal check before each use as specified in the “[Wearing and Adjusting the Mask](#)” section. If you cannot achieve a proper seal, do not use the mask.
- Some filters may lose their filtering ability after gross contamination of the filtering surface with fluids. Please refer to the Instructions for Use of the specific filter that you are using.
- Pneumask PPE can help protect your eyes, nose, mouth, and lungs against certain airborne contaminants. They will not prevent entry through other routes such as the skin, which would require additional personal protective equipment (PPE).
- Pneumask PPE is designed for occupational use by healthcare workers who are

Chat messages:

- Burns, Candy: We don't talk about this anywhere else?
- Abel, Shannon: I think this is a carry-over from when they wanted two version (including a surgical one) – I highlighted as I think it should be removed.
- Gagner, Eric: Definitely remove this. Maybe replace it with a warning to never block the exhale valve. The system will not work well with chin valve blocked.
- Ray Chang: I move it to the "Caution" section and change the wording. I still think this is something that need to be clearly stated in the user manual.
- Gagner, Eric: I don't know the right section but we should have a warning that states something like: discontinue use if the eye chamber becomes fogged or covered in condensation or something. It is normal for the mouth and nose chamber to become fogged.

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https://www.3m.com/3M/en_US/speedglas-welding-helmets-us/adflo/

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“Seal test of the chin valve

A chamber was glued around the chin valve using epoxy, and an 8mm tube was connected to this chamber (**Figure 9-A**). The 8mm tube is connected to an opened water tank (**Figure 9-B**).

First, the water was allowed to equilibrate by manually opening the chin valve briefly (**Figure 9-B**). Then, the mask was lifted up to create a water column of 30mm (**Figure 9-C**), which corresponds to a negative pressure of 30mmH₂O applied to the chin valve. As time passes, the small leakage of the chin valve reduces the negative pressure and therefore the height of the water column.

The leak of the chin valve can be approximated as the ratio between the volume of air in 10mm of tube length, and the time it takes for the water column to go from 30mmH₂O to 20mmH₂O. On average, this time was 21 seconds. As a consequence, in steady state and for a negative pressure of 25mmH₂O on average, the leak of the chin valve is on average 1.5mL/min (12 measurements, standard deviation 0.3mL/min). The leak flow of the exhale valve is lower than the maximum flow allowed to comply with the NIOSH regulation (30mL/min for 25mmH₂O of pressure).


Figure 9. Setup for measuring the leak flow of the chin valve of a Decathlon FreeBreath mask. **(A)** Connection to an 8mm tube. **(B)** Connection to an open water tank - system at equilibrium (chin valve opened). **(C)** The mask is lifted up (chin valve closed) to create a water column of 30mm, which means that a negative pressure of 30mmH₂O is applied to the chin valve.
“

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https://grabcad.com/library/v1-v2-easybreath-decathlon-subea-iso-adapter-rev-7-1/details?folder_id=8121931

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Dear José and Antonio,

hope that you are all doing as you can in this situation

I am sorry but first of all we do not produce PPE we just supply snorkeling mask which is adapted by third party to the mask for the teams in “first line”

All info about this adaptation you can find here 3D plans of adaptator

<https://covmask.cz/en/>

All the drawings of easybreath mask is in the hands of our teams in France as they collect all the adaptations of this mask in the Europe to be able to provide the best possible support and experience .

In copy one of our colleagues Quentin so he might have some answers for you

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<https://antwerpdesignfactory.be/decathlon.html>

A notice on the use of Decathlon Easybreath surface snorkeling masks as Personal Protective Equipment

TLDR; The Decathlon Easybreath snorkeling mask should not be used as a makeshift replacement for FFP2 masks.

PortaCount fit measurements with multiple masks and subjects show that leakage is inherent to the design of the Decathlon mask.

Introduction

The severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) pandemic has led globally to acute shortages of personal protective equipment (PPE) for healthcare workers acknowledged by both governmental and intergovernmental health agencies [CDC, EC, WHO]. A key shortage in PPE for healthcare workers are particulate filtering facepiece respirators, which could either be non-reusable half-face masks, reusable half-face masks or reusable full-face masks. The reusable half- and full-face masks typically use external filter cartridges which are attached to these masks, and these cartridges filter the external air. These filter cartridges are typically classified by their filtration performance, described as the percentage filtration of airborne particles including aerosol droplets that could transport pathogens such as bacteria and viruses, including the COVID-19 virus. In Europe, these cartridges are labelled as P2 (at least 94% filtration) and P3 (at least 99.95% filtration), as per the European Standard EN 143:2000. In the United States, filter cartridges that are not resistant to oil are labelled as N95 (at least 95% filtration), N99 (at least 99% filtration), and N100 (at least 99.97% filtration).

Since about a month a team of product developers has been working on setting up a mass production line for personal protective masks at the University of Antwerp. As everyone knows nowadays the supply of adequately protecting masks are scarce and the health of medical personal cannot be guaranteed without. There are some gradations of masks that are being used in the hospitals: Surgical masks, FFP2 masks and FFP3 masks.

We do not want to make claims about the use or necessity for certain types of masks, but FFP2 is considered the minimum to protect care workers in the vicinity of known sick patients. FFP3 provides a better protection for the care giver but is less comfortable to wear for hours and is equipped with a breath out valve. This return valve is not properly filtered and can pose a risk in the odd situation that the care giver is contaminated, and he can become spreader of the virus.

Surgical masks are a good way of preventing a contaminated wearer from spreading the virus through splashes of biological fluids. These surgical masks do not protect the wearer from becoming contaminated through airborne particles.

As part of the design process, validation was considered a must and thus some onerous basic testing equipment was sought and sourced. But also, an overview of parallel solutions of various other teams were monitored in order to compare and evaluate the advantages and disadvantages. Thus, also the opportunity of the Decathlon Easybreath surface snorkeling mask was checked, and some evaluation models of the mask were purchased days before the end of sales to non-medical instances. As a purely informative way a test was performed to check the fit score of this Decathlon masks compared to a standard FFP2 mask from 3M.

A respirator is only as effective if it fits to the user's face. There are two accepted methods for determining the quality of fit for a user; one is a qualitative test and the other is a quantitative test. In the qualitative tests, the user wears the close-fitting mask and places their head in a large, loose hood into which a scented aerosol is sprayed with a nebulizer. The user performs a series of breathing exercises and head and body movements associated with normal work patterns. If the user can smell the scented aerosol during the exercises, then the face fitting is poor. An example of the qualitative test is the 3M Fit Test Kit FT-30 (bitter) or the FT-10 (sweet), which are typically used by healthcare workers while they are fitting their respirators masks in training. In terms of quantitative tests, there are two relevant testing procedures: the generated aerosol quantitative fit testing protocol and the Ambient aerosol condensation nuclei counter (CNC) quantitative fit testing protocol. In the first protocol, an aerosol is generated from common particulates (for example, sodium chloride) which have sub-micron scaled size distributions. The primary measurement instrument for the second protocol is TSI Incorporated's PortaCount. This device using particulates in ambient air for the testing, but the detection side of the instrument only counts particulates with size distributions between 0.04 – 1.0 μm . In both quantitative fit testing protocols, the users undergo a range of breathing exercises and body movements, from which an overall fit factor can be established. In the United States, the Occupational Safety and Health Administration (OSHA) (29 CFR 1910.134) additionally requires an annual respirator fit test to confirm that the fit of any respirator forms a tight seal on the wearer's face before it is used in the workplace. The test must yield a fit factor that is greater than or equal to 100 for close-fitting (non-)reusable half-face masks, or greater than or equal to 500 for close-fitting full-face masks [OSHA].

We present here, for the first time to our knowledge, quantitative fit factor testing of a full-face snorkeling mask repurposed as emergency PPE for the respiratory protection of healthcare workers. Here, fit factor user testing is performed on a range of Decathlon Easybreath full-face snorkeling masks, 3D printed adaptors and HME filters. Even eliminating the 3D printed adapter from the equation. This letter is separated into the following sections; first a methodology section describes the determination the filtration efficiency of HME filters, the assembly of the full-face snorkeling mask respirator systems and the preparation and testing of the systems with a PortaCount instrument; secondly the results are presented and finally a discussion and conclusion is provided.

Methods

The quality of the Decathlon Easybreath Surface Snorkeling Mask is assessed using a quantitative fit testing method (PortaCount Pro+ 8038, TSI Incorporated). The device measures the number of particles on the outside of the tested respirator as well as on the inside of the respirator and calculates a fit factor (FF) as the ratio between both. In normal use of a face mask, the air that the test person breaths in is filtered and after a few breaths almost none of the particles should be left inside of the mask. If the mask does not seal properly (either on the contact with the face, on any features of the mask, valves, filters or connections), outside air with particles flow straight into the mask and into the sensor. In case of a leak air will not come from the filter that has a relative high pressure drop but will come through the leaks.

Furthermore, an Assigned Protection Factors (APF) is defined, which is the level of protection that should be expected for a respirator. An APF of 10 means that the type of respirator can be safely used in an atmosphere that has a hazardous concentration of up to 10 times the Permissible Exposure Limit (PEL) for that hazard. APF's are determined by the government or a standards organization. In the United States, the National Institute of Occupational Safety and Health (NIOSH) and the American National Standards Institute (ANSI) both establish APF's for various types of respirators. Most full-face negative pressure air purifying respirators (figure 1c) typically have an APF of 50. A half-face negative pressure air purifying respirator (figure 1b) typically has an APF of 10. And Surgical masks (figure 1a) only have an APF of 0.25 to 1. The decathlon Easybreath solutions could be compared to the full-face solutions, but for now will be compared to the average half-face mask standard.

Figure 1 a) Surgical masks, b) A half face negative pressure air purifying respirator, c) full-face negative pressure air purifying respirator.

The accepted practice in the United States is to use ten times the APF as the Fit Factor Pass Level for quantitative fit testing. Considering these APF numbers, the accepted fit factor on a decent mask for medical use in the current COVID-19 crisis can be determined. A Fit Factor Pass Level of 100 (which means a 100 times cleaner air inside the respirator than outside the respirator) and 500 should be used for respectively half and full-face negative pressure air purifying respirators. Any score lower than the pass level means that the mask lets particles enter the mask at such an ease that also viruses can enter.

During the fit testing measurement, the user is instructed to perform several tasks to simulate regular conditions of use of respirators, consisting of seven exercises of one minute: 1) normal breathing, 2) deep breathing, 3) move the head side to side, 4) move the head up and down, 5) talking, 6) bending over and 7) normal breathing again. For each exercise, as well as for the overall measurement, a fit factor is determined.

This study analysed the fit factor for the Decathlon Easybreath mask, which has size S-M (Figure 2b) connected with a heat and moisture exchanger (HME) filter (Altech) (Figure 2c). The filter is connected to the mask using various adapters that are available online (like figure 2a), but in the end is hotmelt glued directly to the mask (Figure 2c) to eliminate doubt about print quality and design variances. Hotmelt glues has also been tested to be a 100% airtight seal for this test method. The quality of this filter was analyzed using the filtration efficiency, which is defined as:

$$Filtration\ efficiency\ (\%) = 100 - \left(\frac{filtered\ count}{ambient\ count} \right) * 100$$

The Decathlon mask is tested in four different forms, to analyze the influence of the adaptor piece and of the valve: 1) without valve(glued shut) (Figure 3a), 2) with normal working valve (Figure 3b), The positioning of the measurement tube was even placed closest to the filter (in the path of the clean air) to be extra certain 3) with valve and repositioned tube (Figure 3c). To provide a comparison, also two masks of 3M are analyzed as well: 1) the Aura 1862+ (FFP2 filter, EN149 classification) and 2) the Aura 1883+ (FFP3 filter, EN149 classification).

Figure 2 a) One of the sample adapters tested, b) One of the tested Decathlon Easybreath masks, b) with accompanying HME filter hotmelt glued to provide a fully airtight seal.


Figure 3 Some of the tested masks; The decathlon Easybreath mask a) with the valve glued shut, b) with the valve still in place, c) with measurement point close to the inlet filter.

Results

To validate the filters quality, the PortaCount measures the particles in the ambient air and compares this to the air coming through the mask. The ambient particles count was 11886 and after the filter 2.44 which results in a filtration efficiency of 99.979%. This filter was then used in all future measurements on the Decathlon masks. Table 1 provides the results of the overall fit factor for the tested respirator models. The difference between the eight test tasks resulted to be negligible for the Decathlon masks and thus only the overall score is shown.

Table 1 The comparison between the Decathlon Easybreath mask against the 3M Aura 1883+ and 1862+ mask. The Decathlon Easybreath is evaluated in three different forms, to investigate the influence of the valve. The minimal required fit factor according to the regulations from OSHA is 500 for full-face respirators as the Decathlon ones but for Covid-19 applications we now use the 100 for half-face respirators, same as the 3M models with FFP2 and FFP3.

Brand	Model	EN standard	Specifications	Overall fit factor / Pass level
-------	-------	-------------	----------------	---------------------------------

Decathlon Easybreath	• With 3D printed adapter	6 / 100
Size S-M	• HME filter	
	• working valve	
Decathlon Easybreath	• Hotmelt glue attached filter	7 / 100
Size S-M	• HME filter	
	• working valve	
Decathlon Easybreath	• Hotmelt glue attached filter	9 / 100
Size S-M	• HME filter	
	• glued shut valve	
Decathlon Easybreath	• Hotmelt glue attached filter	5 / 100
Size S-M	• HME filter	
	• working valve	
	• test point closest to filter	
Decathlon Easybreath	☒ With 3D printed adapter	8 / 100
Size S-M	☒ HME filter	
	☒ glued shut valve	
	☒ With Vaseline added to the silicone to aid seals	
3M	Aura 1862+ EN149	FFP2 filter
		361 / 100
3M	Aura 1883+ EN149	FFP3 filter with valve
		129 / 100

Discussion

The results show that the Decathlon Easybreath does not pass the fit test in the standard form and in all modified forms, with fit factors lower than 10. The 3M masks provide much higher values of the overall fit factor and satisfy the requirements of the Pass level of 100. The magnitude of difference between the Decathlon Easybreath models and the 3M models clearly show the insufficiency of the fit of the Decathlon models. In comparison a surgical mask is required to achieve around 10 on this scale.

During the experiments with the Decathlon mask, the mask was pulled to the face with more than normally advised pressure, up to the point of obvious discomfort for the subject. The silicone part fits well on the face, but the mask is unsuitable for usage in contaminated environments. Furthermore, the difference in fit factor between all modified forms of the Decathlon masks, with and without valve and different tube positioning, is negligible. The filtration efficiency of 99.979% of the HME filter indicates that the inefficiency of this method is to be found in the design of Decathlon mask.

This mask was also tested on different test subjects and on multiple days to eliminate further doubt. But to have clear and comparable results we chose to show the results of a single test subject on all the tests. The tests with different subjects yielded even worse results. So, the data we show here is the best we could get out of the mask.

These measurements give us reference values of the commercially available Decathlon mask. Optimizations regarding the filtration and fitting were applied to analyze the effect thereof. However, the focus of this study is not to design or modify the commercial model, since we are focusing on designing qualitative FFP2 and FFP3 masks that we can produce in large quantities and use the obtained results to compare with our design solution.

We as a team are all great fans of 3D printing and modifying tools for other purposes, but fear that the use of these mask will provide a false feeling of protection to anyone using them. In our tests the decathlon mask should be looked at as an alternative for surgical masks rather than for personal protection equipment like the FFP2 and FFP3 masks provide.

Designing other adaptor pieces, using different materials or different 3D print settings will not aid in improving the results that we measured since we took the 3D printed part out of the equation.

Only the switch to a positive pressure filtered mask design could be a possibility. This way there is a constant airflow from the filter's through the mask. Leaks would in that case still not allow contaminated air to creep into the mask. But only at sufficient airflows is this certain and thus also here we only advise to continue this path if minimal testing can be done to verify the safety aspects.

Conclusion

The fit test shows that the Decathlon mask is not an airtight mask and thus should not be used as a solution for filtration as an alternative to FFP2 and FFP3 masks. Based on our tests, the issue is not in the filter cannister and not in the 3D printed adaptors. The leakage is inherent to the design of the Decathlon mask. The mask was designed for a different application where the outside pressure of the water helps to seal the various parts of the design to each other. And this mask will require intensive redesigns to achieve a decent airtight fit.

As a design team we did not plan on doing the extensive set of tests on the Decathlon mask that we did in the past couple of days. We just had purchased a couple of the masks in a fling in case this would become the solution for all care workers. But as soon as we had the test, we felt the need to see how it would compare to the standard and to our own design. Expecting an off the chart result given the look of a full-face mask we were shocked to see so low of a score. Even lower than pretty much any surgical mask would receive. We took apart some of the masks, tightened all the joints, even used Vaseline all around the perimeter to help the mask provide a decent seal. Thus, our advice is to use this mask for nothing else than you would feel safe to do with a simple surgical mask and a face shield as this solution does not provide security to airborne particles like an FFP2 mask should.

References

[CDC] National Center for Immunization and Respiratory Diseases (NCIRD), Division of Viral Diseases “Strategies to Optimize the Supply of PPE and Equipment” Retrieved from <https://www.cdc.gov/coronavirus/2019-ncov/hcp/ppe-strategy/> (published 3 April 2020, last accessed 4 April 2020)

[EC] Balazs UJVARI “COVID-19: Remarks of Commissioner for Crisis Management Janez Lenarčič at the midday press briefing on the rescEU stockpile” Retrieved from https://ec.europa.eu/commission/presscorner/detail/en/STATEMENT_20_495 (published 19 March 2020, last accessed 4 April 2020)

[WHO] Fadela Chaib, “Shortage of personal protective equipment endangering health workers worldwide” Retrieved from <https://www.who.int/news-room/detail/03-03-2020-shortage-of-personal-protective-equipment-endangering-health-workers-worldwide> (published 3 March 2020, last accessed 4 April 2020)

[OSHA] Occupational Safety and Health Administration (1998). Respiratory protection (Standard No. 1910.134). Retrieved from <https://www.osha.gov/laws-regs/regulations/standardnumber/1910/1910.134> (last accessed 4 April 2020)

F.A.Q.

Question

Answer

But my design of the valve is different in this or that way:

We have tested multiple adaptor pieces and found the problem not to lie in the adaptor pieces but rather in the design of the mask itself. Thus, no adaptor (even injection molded) could improve this seal.

I have a different mask/solution for decent FFP2 masks, can I use your testing setup?

We can't offer that service I'm afraid. We are a design team working on a solution of our own on a big scale. We sourced the device to do our in-house tests so that we can run dozens of tests a day on our design iterations and providing such a service for everyone is just impossible. The results we noticed on these Decathlon masks were just so obviously bad that we felt like we had to share this as soon as possible. But to make sure we did not miss any obvious fault we spend more than a day trying alternative methods and modifications. Without luck.

You probably just had a bad example of the mask!

We tested 5 different masks in total. The ones you see in the pictures are in fact various examples of the same mask. Just to be certain and to evade this possibility of a bad example. We also tested both sizes S-M and M-L.

There is a new version of the Easybreath mask! Did you test that one as well?

The masks we acquired were purchased on the 10th of March but are of the old design. We have not had our hand on the new design, but from what we gather it is not fundamentally changed in a way that would affect our measurements.

Is this report also on the use of Decathlon masks for positive pressure systems?

We have not done any tests of these masks with a positive air pressure system for two reasons; First because according to the regulations a positive air system should be tested with the pressure system disabled. The PortaCount device pumps air inside by itself. On the other hand, positive air pressure systems could work since all air inside the mask is flushed with clean air constantly. But when the user breathes heavily we are unsure that the ventilator can keep up at a sufficient pace. Thus, we recommend anyone working on such solutions to investigate the regulations and to at least do some basic validation with quantitative measurement techniques.

We are also working with these masks and have done fit tests as well; we have way better results!

We would love to hear from you with a little report of the used material and methods. If this is so then we are more than willing to update or even take down this report. But given the measurements we have at the moment we feel like no proper analysis had been conducted on these systems yet.

How can we get in touch?

Our mailbox and phones are exploding nowadays with questions, therefore we provide this F.A.Q. overview here above. If you feel positive you have a good update for us or prove of better tests on this mask then you can contact us at: decathlon@antwerpdesignfactory.be.

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Cleaning Instructions

****The following instructions are based on the instruction of OSHA on cleaning a respirator. Please consult EH&S officers or the Infection Control Center of your hospital to achieve a consensus of appropriate cleaning SOPs.****

1. Perform hand hygiene.
2. Put on a gown, gloves, and a protective mask.
3. Use alcohol-based wipe or hydrogen peroxide wipe to wipe the surface of the mask and filter. Gently disconnect the filter from the mask. Dispose the wipe as biohazardous waste.
4. Dispose the filter if it is designed for one-time usage, or put it into the dedicated area if the filter is designed for disinfection and reuse.
5. Wash mask and filter adapter thoroughly in warm water with a mild detergent.
6. Rinse the mask and adapter thoroughly in clean, warm, running water. Drain.
7. Prepare a 50ppm chlorine solution: add one milliliter of laundry bleach per liter of warm water. Prepare enough solution to fully immerse the mask.
8. Immerse the mask and adapter in a hypochlorite solution (50 ppm of chlorine) for 2 minutes.
9. Wearing fresh gloves, rinse the mask and adapter thoroughly in clean, warm, running water. The mask and adapter must be thoroughly rinsed with water to remove any detergents or disinfectants that may result in dermatitis or damage of the mask components.
10. Dry the mask and adapter with a clean lint-free cloth or allow to air dry.
11. Reassemble the mask and store in a clean space.

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https://1drv.ms/v/s!AshugVzknI_vh74WRcYPm6tE59e-gA?e=HQqFYc

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


		etCO	FiCO	F	Nota
Peça em Y twisted	2		2	R	
Válvula Aberta para o exterior. Um filtro em cada ramo	22:5			1	Mask is comfortable. No changes to the original design. Little fogging.
	5	42	10	2	
Sensor no interior da máscara	23:0			1	Comfortable. Little fog. Air flows in both entries. Adjusted the mask. The seal around the nose is not perfect.
	0	43	7	3	
	23:0			1	
	5	42	10	4	
	23:1			1	Itching in my nose. Not much fog
	0	39	7	5	
	23:1			1	
	5	40	5	5	
	23:2			1	
	0	40	8	5	
	23:2			1	I want scratch my face badly. Not much fog
	5	40	10	6	
	23:3			1	
	0	42	9	4	
	23:3				
	5				
	23:4				
	0				

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Registos dos testes de CO2

https://1drv.ms/x/s!AshugVzknI_vh75B9XvV3OAmqBqX2g?e=VX0IQg

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
Referências

Massimiliano Sorbello, MD
Anesthesia and Intensive Care
AOU Policlinico Vittorio Emanuele
Catania - Italy
Postgraduate School Anesthesia and Intensive Care University of Catania
SIAARTI Airway Management Study Group Task Force - Coordinator
European Airway Management Society Board Member
(H) +39 (0)95 7435268
(M) +39 349 6277107

Outros links de interesse

<https://www.materialise.com/en/manufacturing/materials/pa-12-sls>

https://www.materialise.com/system/files/resources/materialise_datasheets_def_file_PA-AF.PDF

<https://grabcad.com/filip.kober-1/models?page=1>

https://grabcad.com/library?page=1&time=all_time&sort=recent&query=snorkel

https://grabcad.com/library?page=1&time=all_time&sort=recent&query=easybreath

<https://forums.autodesk.com/t5/fusion-360-design-validate/y-piece/td-p/9427916>

<https://www.thingiverse.com/tag:covid-19>

<https://www.thingiverse.com/search?q=snorkel+mask&type=things&sort=relevant>

https://drive.google.com/drive/folders/14Q3TEl5JVeN2QpDpKo1Alx_wnGeolKlK

Um miniventilador DYI

<https://www.instructables.com/id/COVID-Swim-Mask-HEPA-PAPR/>